

Congo Red Stain Kit

PRODUCT INFORMATION:

REF

SSP002 100ml
SSP002 250ml
SSP002 500ml

PERFORMANCE CHARACTERISTICS:

Staining Interpretation:

Nuclei: Dark Blue

Amyloid:

In Transmitted light – Pink to Red

In Polarized light – Apple green birefringence

Connective tissue; Collagen: Light Red

SUMMARY AND EXPLANATION

For laboratory use only

The Congo Red Stain Kit is intended for use as a qualitative histological stain to demonstrate Amyloid in formalin-fixed, paraffin-embedded tissue selectively. This product is not intended for diagnostic or therapeutic use. The results are to be interpreted by qualified personnel in conjunction with other clinical and laboratory findings.

PRINCIPLE OF THE PROCEDURE

Amyloid is a homogenous structure made up of protein fibrils (each between 8 and 15nm in diameter) that can be stained eosinophilically, e.g., in the case of amyloidosis. Amyloid deposits are formed in the intercellular spaces. All deposits of amyloid contain similar protein fibrils that are resistant to the body's natural defence mechanisms, and once they have formed, cannot be eliminated.

Congo Red is a metachromatic anionic dye, and it forms the hydrogen bridge bonds with the carbohydrate component of the substrate. The beta pleated sheets of amyloid are suitable in size and shape to accommodate the Congo red molecules, which are held in the lattice work of the beta pleated sheets. The tissue stained with Congo red appears pink-red under the transmitted-light microscope and a brilliant apple-green birefringence under the polarized light due to conspicuous dichroism. Congo red may also stain unexpected structures such as keratin, elastic and dense collagen fibers, which cannot be visualized under polarized light.

REAGENTS PROVIDED

Kit Contents	Product Code	Storage Conditions	Pack sizes		
			100ml	250ml	500ml
Congo Red Stock Solution (Reagent A)	IPS019	RT	100ml	250ml	500ml
NaOH - A (Reagent B)	IPS024	RT	2ml	5ml	10ml
Alkaline Alcohol Solution (Reagent C)	IPS020	RT	100ml	250ml	500ml
Modified Mayer's Hematoxylin (Reagent D)	PS020	RT	100ml	250ml	500ml

STORAGE AND HANDLING

Storage Recommendations: Store at Room Temperature. When stored at the appropriate conditions, the reagents are stable until expiry. **Do not use the reagents after the expiration date provided on the vial.**

To ensure proper reagent delivery and stability, replace the dispenser cap after each use and immediately place the bottles at room temperature, away from sunlight, in an upright position.

During transport, short-term exposure to 2- 8 °C does not affect product performance.

PREPARATION OF WORKING SOLUTION

For 10ml of Congo Red Working Solution: Measure 9.9ml of Reagent A (Congo Red Stock Solution) and add 0.1ml of Reagent B (NaOH - A).

Laboratory Use Only

Note: Filter the solution and use it 20 minutes before it may degrade in prolonged standing.

SPECIMEN PREPARATION

Recommended positive controls:

Formalin-fixed, paraffin-embedded or frozen tissue sections containing amyloid deposits.

Sample preparation and fixation:

- Fixation plays an essential role in preserving the tissue structure to be visualized using the stain.
- 10% Neutral buffered formalin can be used for fixation.
- Ensure that the fixed sections are adequately embedded in paraffin.
- Cut the sections, usually 6-12 µm, to show smaller amyloid deposits which can exhibit birefringence under polarized light.

PRECAUTIONS

- Normal precautions carried out in handling laboratory reagents should be followed
- This product should be used by qualified and trained professional users only
- The product contains alcohol and is classified as highly flammable; it must be kept away from ignition sources
- It can cause serious eye and skin irritation. Refer to Material Safety Datasheet for any updated risk, hazard or safety information
- Dispose of waste, observing all local, state, provincial or national regulations
- Do not use reagents after the expiration date
- Use protective clothing and gloves while handling reagents
- Avoid contamination of reagents, as it may lead to incorrect results

MATERIALS REQUIRED, BUT NOT PROVIDED

- Xylene
- Graded Alcohols (50%, 70%, 95%, Absolute)
- Bleuing Solution (Optional)
- DPX Mountant
- Microscopic Slides (Positively charged)
- Slide Holder
- Polarized Lens
- Cover Slips
- Coplin Jars

STAINING PROCEDURE

- Bake the FFPE tissue sections at 70°C in a hot air oven for 20 minutes.
- Deparaffinize and hydrate the sections using xylenes, graded alcohols - 80%, 70%, 50% and distilled water for 2 minutes each.
- Stain the sections using Congo Red working solution (Refer to preparation of working solution) for 10 minutes.
- Rinse in distilled water for 1-2 minutes.
- Differentiate quickly in Alkaline Alcohol solution (Reagent C) 9-10 seconds.
- Rinse in tap water for 2 minutes.
- Counterstain with Modified Mayer's Hematoxylin (Reagent D) for 30-60 seconds.
- Rinse the slides under running tap water for 2 minutes.
- Place the slides in distilled water for 10 minutes.
- Dehydrate rapidly using graded alcohols.
- Clear the sections using two changes of xylene.
- Cover slips the sections using DPX mountant.

QUALITY CONTROL

The recommended positive tissue control for the Congo Red Stain Kit is tissue sections containing amyloid deposits.

PERFORMANCE CHARACTERISTICS

Congo Red Stain Kit stain for Nuclei in blue colour, Amyloid - In transmitted light appears Pink to Red in colour and in polarized light appears Apple green birefringence in colour and Connective tissue and Collagen stains Light Red in colour.

TROUBLESHOOTING

1. Follow the specific protocol recommendations according to the data sheet provided
2. Tissue staining is dependent on the handling and processing of the tissue prior to staining. Improper fixation, tissue processing, freezing, thawing, washing, drying, heating, sectioning or contamination with other tissues or fluids may produce artifacts, reagent trapping or inaccurate results
3. Do not allow the section to dry out during the entire staining process
4. Gently mix all the reagents prior to use.
5. Excessive or incomplete counterstaining may compromise the interpretation of the results
6. If unusual results occur, contact PathnSitu Technical Support at +91-40-2701 5544 or E-mail: techsupport@pathnsitu.com

LIMITATIONS AND WARRANTY

1. This product is intended for use only by authorised, trained, and qualified personnel.
2. A qualified and trained pathologist/personnel must interpret the results of the test.
3. Interpretation of test results must be made in conjunction with relevant background information and additional laboratory findings.
4. Always use the recommended volume and concentration of reagents to ensure complete coverage of the tissue section and to minimise the risk of false-positive or false-negative results.
5. Use appropriate buffers, instruments, consumables, and incubation conditions as recommended to achieve optimal staining performance.
6. It is strongly recommended to include known positive and negative controls when performing the test to ensure the validity of results.
7. The product has been validated on formalin-fixed, paraffin-embedded (FFPE) tissues. The end user must establish performance on other tissue types.
8. Unexpected results may occur in untested tissues due to inherent variability in tissue components.
9. False-positive reactions may occur due to insufficient washing, inappropriate protocol conditions, or other contributing factors.
10. In instances where the staining pattern or localisation differs from the specifications outlined in this datasheet, please get in touch with technical support for guidance.
11. Maintain the product under the recommended storage conditions to preserve reagent stability and performance.
12. Do not use reagents that appear cloudy, discoloured, or show signs of contamination. Discard any components showing signs of deterioration.
13. NaOH is strongly caustic and unstable in air (absorbs CO₂); incorrect concentration alters Congo Red binding and reduces staining quality.
14. This product is intended for single-use application only. Once applied to a tissue section, reagents should not be recovered or reused, as this may compromise test integrity and specificity.
15. PathnSitu makes no warranties beyond those expressly stated in the product description.
16. PathnSitu shall not be liable for property damage, personal injury, time or effort, or economic loss arising from the use of this product.
17. Please refer to the complete datasheet for all instructions, precautions, and additional product limitations.
18. For detailed information and specifications on individual components, please refer to Product Material Safety Data Sheet (MSDS)

BIBLIOGRAPHY

1. Romeis – Mikroskopische Technik, Editors: Mulisch, Maria, Welsch, Ulrich, 2015, Springer-Verlag Berlin Heidelberg.
2. Theory and Practice of Histological Techniques, John D Bancroft and Marilyn Gamble, 6th Edition.
3. Selected Histochemical and Histopathological Methods. Edited by SW Thompson Charles C, Thomas, Springfield (IL), 1966, pp 402-405.
4. Theory and practice of Histotechnology. Edited by DC Sheehan and BB Hrapchak, 2nd ed. Mosby. St.Louis (MO), 1980, pp177-178.
5. Conn's Biological Stains: A Handbook of Dyes, Stains and Fluorochromes for use in Biology and Medicine, 10th Edition. (ed. Horobin, R.W. and Kiernan J.A), Bios, 2002.

EXPLANATION OF SYMBOLS

LOT



Lot number / Batch number

Storage limitation

Date of manufacture

Manufacturer address



Expiry

RT

Room Temperature

REF

Catalogue number