

# TISSUE PROCESSOR TP3S1



## Features

- Capacity of 200+ cassettes holding tissue basket capacity
- Dual Robotic Arm for increased throughput
- Equipped with self-rechargeable battery
- Bright LCD display for easy visibility
- User can set multiprogramming capacity
- User friendly software
- Robust instrument

# Description

Tissue processing comprises of majorly 3 steps (Dehydration, Clearing & Impregnation) using multiple reagents which is required to take an animal/human tissue from fixation to the state where it is completely infiltrated in an appropriate histological wax and can be readily embedded for section cutting using a microtome.

## Technical Specifications

### 1. Power Voltage: AC 220V-240V

Frequency: 50±10%Hz

Power: 550W

### 2. Basic Parameters and Performance

Hydraulic Jar quantity	14 units (10 Chemical Jars and 4 Paraffin Jars)
Single Jar Volume	2000ml
Single Jar Treatment time	Jar 1 can be freely set within 0~99 hours. The other Jar could be freely set within 0~24 hours
Temperature Presetting Range for Wax Jars	Room temperature to 80°C
Residual Liquid dropping time while replacing Jar	Approx 30 seconds
Dimensions (in mm)	1370(L)×440(W)×525(H)
Total Weight	Approx 102 kg

## Ordering Information

Equipment Name	Catalog No
Tissue Processor	PS-TP3S1

## Contact Us

### USA

1257, Quarry Ln, Suite150  
Pleasanton, CA-94566, USA  
Phone: +1 925-218-6939

### Email

techsupport@pathnsitu.com  
customerservice@pathnsitu.com

[www.pathnsitu.com](http://www.pathnsitu.com)

### India

PS Heights, 2nd Floor, Plot# A-28/1/11B/B,  
Road#15 Nacharam IDA,  
Hyderabad-500 076, Telangana.

Phone: 040-2701-5533 /5544  
Fax: 040-2701-5544

