

Microscopic Slides-Colored

PRODUCT INFORMATION:

PS7109-50

50 slides per pack

INTENDED USE

The Microscope Slides are colored on one end to help ensure specimen identification and traceability. The colored ends are resistant to stains and chemicals commonly used in laboratories. The microscopic slides are all pre-cleaned and ready to use which may be used across a wide range of applications which include Pathology, Cytology, Microbiology and Hematology laboratories which accommodates various marking methods such as conventional labeling, writing, or printing. With multiple colors, the microscopic slides are thoughtfully tailored to suit the laboratory needs.

SUMMARY AND EXPLANATION

The Coloured Microscope slides have one edge surface with a pink, blue and green label field. Hence different preparations can be clearly differentiated by their label. The color printing of the label field is resistant against all solvents normally used in the laboratory. These are pre-cleaned clear glass slides with grounded edge of 25.4 x 76.2mm (1" x 3") and (1mm - 1.2mm) thickness which help maintain specimen integrity and achieve accurate results.

STORAGE AND HANDLING

The coloured microscopic slides should always be stored in a dry place, to avoid sticking and corrosion of the glass surface. The changes in temperature and humidity results in moisture contamination of the slides. Slides should not be stored on concrete floors and should be kept as far away as possible from heat and air conditioning to avoid large variations in temperature and humidity. Do not store slides close to solvents as this might turn their surfaces hydrophobic.

RECOMMENDED PROTOCOL

- Slides are normally ready for use.
- Labelling fields on the slides serve to clearly identify the preparations, so they must be labelled appropriately.
- The slides can be used according to the respective laboratory validated procedures.
- Avoid mishandling of slides to minimize scratches and breakage.
- PathnSitu coloured microscopic slides are developed to work with pathology, cytology, microbiology and hematology laboratories of any source.

PRECAUTIONS

- The slides are intended for single use by qualified personnel only.
- They must not be used more than once and, after use they must be properly stored/disposed.
- Wear gloves to reduce the risk of injury from sharp glass edges as much as possible
- Any broken slides must be discarded into containers which are purpose-built for sharps.

TROUBLESHOOTING

Follow the reagent specific protocol recommendations according to data sheet provided. If unusual results occur, contact PathnSitu Technical Support at 040-2701 5544 or techsupport@pathnsitu.com

LIMITATIONS AND WARRANTY

There are no warranties, expressed or implied, which extend beyond this description. PathnSitu is not liable for property damage, personal injury, or economic loss caused by this product.

EXPLANATION OF SYMBOLS

CAT# - Catalogue number

Qty: Quantity of box / No. of slides per box