

### Microscopic Slides- Plain slides

#### PRODUCT INFORMATION:

PS7101-72

72 slides per pack

#### INTENDED USE

The Plain microscopic Slides are high quality pre-cleaned glass slides which are ideal for use in all routine laboratory applications. The plain microscopic slides are all ready to use which may be used across a wide range of applications which include Pathology, Cytology, Microbiology and Hematology laboratories which accommodates various marking methods such as conventional labeling or writing.

#### SUMMARY AND EXPLANATION

Plain Microscope Slides provide a clear, classic and even surface produced from the finest quality of white glass. The plain slides are per cleaned clear glass slides with grounded edge of 25.4 x 76.2mm (1" x 3") and (1mm - 1.2mm) thickness which help maintain specimen integrity and achieve accurate results.

#### STORAGE AND HANDLING

The plain microscopic slides should be always stored in a dry place, to avoid sticking and corrosion of the glass surface. The changes in temperature and humidity results in moisture contamination of the slides. Slides should not be stored on concrete floors and should be kept as far away as possible from heat and air conditioning to avoid large variations in temperature and humidity. Do not store slides close to solvents as this might turn their surfaces hydrophobic.

#### RECOMMENDED PROTOCOL

- Plain Slides are normally ready for use.
- Labelling the slides serve to clearly identify the preparations, so they must be labelled appropriately as the plain slides are designed with no labeling field.
- The slides can be used according to the respective laboratory validated procedures.
- Avoid mishandling of slides to minimize scratches and breakage.
- PathnSitu microscopic plain slides are developed to work with pathology, cytology, microbiology and hematology laboratories of any source.

#### PRECAUTIONS

- The slides are intended for single use by qualified personnel only.
- They must not be used more than once and, after use they must be properly stored/disposed.
- Wear gloves to reduce the risk of injury from sharp glass edges as much as possible
- Any broken slides must be discarded into containers which are purpose-built for sharps.

#### TROUBLESHOOTING

Follow the reagent specific protocol recommendations according to data sheet provided. If unusual results occur, contact PathnSitu Technical Support at 040-2701 5544 or [techsupport@pathnsitu.com](mailto:techsupport@pathnsitu.com)

#### LIMITATIONS AND WARRANTY

There are no warranties, expressed or implied, which extend beyond this description. PathnSitu is not liable for property damage, personal injury, or economic loss caused by this product.

### EXPLANATION OF SYMBOLS

**CAT#** - Catalogue number

**Qty:** Quantity per box / No. of slides per box