

# FLEXIGROSS

## SMART GENERATION GROSSING



FOR DETAILS SCAN



**ENHANCED USER SAFETY**

**MACRO DIGITAL DOCUMENTATION**

**ADJUSTABLE HEIGHT**

**TAILOR-MADE MODULARITY**

**OPTIMAL ERGONOMICS**

**ADJUSTABLE WORK AREA**

PathnSitu's flexigross series is specially designed for histology labs, providing sturdy, customizable, and affordable workstations. flexigross series workstation supports for higher suction and additional features to assist user in grossing, while also protecting user from harmful formalin fumes.

At PathnSitu, we understand that every lab has unique requirements, which is why our flexigross series can be tailor-made to fit to user specific needs. The series features medical grade stainless steel which slide across the entire length of the basin underneath. Unique feature of this instrument is that it comes with the technology to adjust the height of the instrument. Along with this it comes with additional unique features such as adjustable work station and wheels for easy mobility. This innovative design can accommodate accessories like Formalin monitor, Label Printer, Camera, Magnifying Glass etc... providing the flexibility to user to customize the workstation.

Our commitment to quality and customer satisfaction is paramount, and we strive to provide reliable and efficient products. PathnSitu's flexigross series for a dependable, customizable, and affordable solution for your grossing needs.

# SPECIFICATIONS

1. Designed for grossing surgical /histology specimens
2. Inbuilt blower with adjustable speeds, capable of sucking out the vapors
3. Grossing station body is built with medical grade stainless steel for durability and easy cleaning
4. Adjustable height from 850 to 1250mm for user comfort and convenience
5. Availability of dual taps and faucets supports for Hot & Cold with foot operations, adjustable hose and 2 different tap flow patterns
6. Adjustable work area from left to right with high-performance sink
7. Larger work area as per the grossing station size for user convenience
8. Works with 220V single phase AC supply for efficient operation
9. Magnetic toolbar to attach grossing tools
10. High-performance LED lighting for better visibility
11. Availability of foot switch operated formalin pump & water taps
12. Mobility: Easy movable with castor wheels for easy relocation



# OPTIONAL ACCESSORIES

1. Camera
2. Computer
3. Label Printer
4. Magnifying Glass
5. Mic Set
6. Soap Dispenser
7. Monitor stand holder
8. End-to-End Histology Solution Software

# ORDERING INFORMATION

Catalog No	Work Area	Catalog No	Work Area
PFGS-3	3 Feet	PFGS-7	7 Feet
PFGS-4	4 Feet	PFGS-8	8 Feet
PFGS-5	5 Feet	PFGS-9	9 Feet
PFGS-6	6 Feet		

Adjustable Height Range	
Minimum Working Height	850 mm
Maximum Working Height	1250 mm



**CONTACT US TO MAKE YOUR OWN  
CUSTOMIZED AND TAILOR MADE  
GROSSING STATION**



## Contact Us

- 📍 **USA** : 7026 Koll Center Parkway, STE#211, Pleasanton, CA 94566, USA. ☎ +1 925-218-6939
- 📍 **Canada** : 30 Bertrand Ave, Unit B209, Scarborough , Ontario - Zip - M1L 2P5. ☎ +1 437-269-3912
- 📍 **India** : PS Heights, Road#15 Nacharam IDA, Hyderabad-500 076, Telangana. ☎ +040-2701-5533 /5544
- ✉ info@pathnsitu.com | customerservice@pathnsitu.com
- 🌐 [www.pathnsitu.com](http://www.pathnsitu.com)